

# Bottle Rocket Launcher

#17000

SB21429m

## OPERATING INSTRUCTIONS

### Purpose:

To launch rockets using commonly available materials such as plastic soda bottles and cardboard tubes. It also provides an exciting introduction to aerodynamics and rocketry.

### Contents:

- One (1) Launcher (includes, locking mechanism and check valve)
- One (1) Plumb line with weight.

### Required Materials:

- Bicycle or other air pump.
- Tire pressure gage
- Dowels or straight sticks (3 or 4)
- Tape
- Glue
- Scrap cardboard
- Water
- Common pin or nail

## Warning!

### Do:

- Do Wear Protective Eye Gear
- Do Wear Protective Clothing (rain gear)
- Do Stake Launcher Securely into the Ground

### Do Not:

- Do Not Aim at Anyone
- Do Not Pressurize Beyond 40 psi
- Do Not Stand Over Pressurized Rocket
- Do Not Stand Within 10 Feet of Launching Rocket

### Assembly:

Congratulations! Your new bottle rocket launcher will let you use any standard plastic soda bottle with capacities of 500 ml to 2 liters! When choosing a bottle to use for rocketry, remove the cap and check the bottle's neck ring for any cracks or defects that may fail when pressurizing the bottle. If the bottle looks good you're ready to begin building your rocket. You're probably anxious to get started right away; and there's no reason you can't! Your bottle should be ready to use, but first remove the older bottle's flat plastic base if it has one. This base is lightly glued onto the bottle and can be pulled off with a few good tugs. Partially fill the bottle with water, according

to the details below. Press the rubber stopper firmly into the neck of the bottle. Notice that each arm of the launcher has a slot cut into its side. Swing the arms up so that each slot fits over the neck ring on your bottle. When the neck ring is firmly held in each slot in each arm, insert the locking bar. See the assembly diagram for details. When you've found an appropriate launch site, press the two long pointed legs of the launcher a few inches into the ground. This will anchor the launcher in place. A quick launch pad can be made by pushing three dowels or sticks into the ground equally spaced around the bottle. These dowels will help hold the bottle upright while pressurizing and launching.

*Nasco*

Fort Atkinson, Wisconsin  
FREE PHONE ORDER SERVICE  
Toll-Free 1-800-558-9595  
Fax 920-563-8296  
www.eNASCO.com  
E-mail: info@eNASCO.com


Before you begin to pump up the bottle, there are a few precautions you must take, and most of these are common sense.

- 1) While pumping up the bottle, periodically check the pressure. Use a pump with a built in pressure gage and measure the pressure while pumping.
- 2) Do not exceed 40 psi.
- 3) Don't look down on the rocket during or after pressurization.
- 4) Don't aim your rocket at any people, buildings, or other structures.
- 5) Always set up your rocket well away from any trees, telephone wires, or houses.
- 6) Do not launch your rocket in parking lots where it may damage someone's car.

Remember, when you pull the cord, you are the pilot in command and are responsible for that particular flight! So be careful! Now that you understand the precautions that must be observed with every flight, the bottle is now ready to be pressurized!

Begin by fueling your rocket with ordinary water. This works because every bit of water that gets forced out of the bottle produces an equal and opposite force that pushes the bottle upward! This "thrust" is proportional to the speed with which the water escapes the bottle.


When the bottle is half full, attach your bicycle pump to the valve stem into the end of the pressurizing hose and start pumping. When you have reached 40 psi you are ready to launch. Stand away from the rocket and tension the trigger string. After announcing the countdown to launch (every professional rocket pilot should count down from 10...9...8...) give the trigger line a sharp tug. This will release the locking mechanism and the rocket will blast off from the pad with lightning speed! This marks the successful launch of your first rocket!



**Figure 1**  
Attaching the rocket launcher to your bottle

#### Discussion:

Now that you've launched some rockets, your head may be full of ideas for better designs to get it to go higher, fly further, or have increased thrust. The following are a few suggestions to help you get started.



**Figure 2**  
Building a launch pad.

The stability of the rocket depends upon the rocket's length. You'll find that longer rockets fly straighter than shorter rockets. Long cardboard tubes make excellent bodies for bottle rockets. The larger diameter tubes often found in rolls of gift and Christmas wrapping paper make the best bodies. Smaller diameter tubes will work but they are not as strong and may not hold up as well under repeated use. Another method of constructing rocket bodies is to tape several paper towel tubes together.

Fins can also improve the stability of the rocket. Experiment with different size and shape fins. The fins do not have to be very large to be very effective! The fins can be made using scrap cardboard from old boxes. Select your cardboard carefully, making sure that it is reasonably flat and does not have bends or folds in it. Draw your fin pattern on a sheet of paper and then use this pattern to cut the number of fins you need from the cardboard sheets. These fins can then be glued or taped to the body of the rocket. The fins could also be taped to the bottle itself! The fins work the best when placed near the bottom of the rocket, but you can also experiment by placing the fins at varying positions along the body.

The nose of the rocket is another point of consideration. The nose of the rocket must do two things: first it must streamline the rocket so the air will pass over it more efficiently, second: it must hold the parachute or other recovery system or it must be able to absorb the impact of the rocket as it hits the ground. The nose cone can be made of cardboard rolled into a cone and then taped to the top of the rocket, or you may want to use a smaller plastic bottle turned upside down so the smooth round bottom of the bottle points upward on the rocket.

The parachute can be made by cutting a circular patch from a large plastic garbage bag. Six lengths of string each about sixty centimeters long (2 feet) can be cut and taped around the perimeter of the parachute. After each string is taped in place, gather the free ends of the string together and tie them into a knot. This knot can then be taped to the inside of the rocket body tube. The parachute can then be gently rolled up and placed in the rocket body. After the parachute is loaded the nose cone can be rested on top of the rocket. Don't tape the nose cone to the rocket body if you want the parachute to come out!

### **Rocket Safety:**

Now that you have finished building your high performance rocket, it is worth repeating a few of the safety rules that have already been mentioned.


- 1) **Do not exceed 40 psi while pressurizing the bottle.**
- 2) **Don't look down on the rocket during or after it is pressurization.**
- 3) **Don't aim your rocket at any people, buildings, or other structures.**
- 4) **Always set up your rocket well away from any trees, telephone wires, or houses.**
- 5) **Do not launch your rocket in parking lots where it may damage someone's car.**

**Remember - as pilot in command, you are responsible for each and every flight!**

### **Rocket Experiments:**

"Whoa! How high did that thing go?" has probably been a question you've asked several times as you watched your rocket soar overhead. There is a relatively simple method you can use to measure the altitude of your rocket! It is simply a sheet of graph paper mounted on a sheet of cardboard with a plumb line attached to one corner.

To assemble the altimeter, paste the altimeter face (found at the end of this manual) to a flat sheet of cardboard that has been cut to the same size. Make sure the top edge of the cardboard is flush with the top edge of the altimeter page. Use a common pin to poke a hole through the cardboard at the top front corner where indicated. Slip one end of the plumb line (thread) through this hole.



**Figure 4**  
Assembling the altimeter.

Tape the thread to the cardboard backing where it comes through the hole. Tie weight (large washer) to the free end of the thread as in Figure 4. Your altimeter is now ready to use.

Let's start by measuring the height of something you know. First, measure the distance between the object you're measuring and the position you will be standing when you measure its height. This distance is called the "BASE" and should be at least 10 meters. Hold the altimeter card at eye level and sight along its top edge. The place where the thread is taped to the cardboard should be away from you. Tip the cardboard upward until the top edge of the cardboard is pointing at the top of the object you're measuring. When the string is hanging freely and no longer swinging, pinch it against the face of the altimeter with your finger, as in Figure 5. You can now read the altimeter by finding an intersection on the graph paper that the string crosses. Read the quantity "Y" from the horizontal axis of the graph paper (this is proportional to the height of the object) and read the quantity "X" from the vertical axis of the graph paper (this is proportional to the distance between you and the object). Using similar triangles, the rocket altitude or "HEIGHT" can be calculated by using the following formula:

$$\text{HEIGHT} = \text{BASE} * (\text{Y}/\text{X})$$

Where:

HEIGHT = height of the object

BASE = horizontal distance from the object

Y = measured from the altimeter (horizontally)

X = measured from the altimeter (vertically)


Don't forget to add your own height to the calculated height! Actually, you should measure the distance from the string hole down to the ground, and add that

to your calculated height. Your most accurate measurements of height will be when your horizontal distance from the object is close to the actual altitude of the object.

The altitude of the object will have the same units as the horizontal distance from the object. In other words, if you measured the distance from the launch pad in meters, then your calculations will produce an altitude in meters. If you measured the distance from the launch pad in feet, then the altitude will be in feet.

Now that you are able to get accurate measurements of the rocket's altitude, you're ready to do some serious experimenting! Remember - change only one thing at a time; keep everything else constant!

For example, to measure the effect of pressure on the rocket's altitude, use the same rocket, pour the same amount of water into the bottle, and then, for each experiment, pressurize the bottle to different amounts. Record each launch on a data sheet that might look like Table 1.



**Figure 5**  
Reading the altimeter

Table 1.

<b>PRESSURE TEST:</b>				
<b>Launc h #</b>	<b>Ounces of Water Used</b>	<b>Rocket Weight</b>	<b>Pressure Used</b>	<b>Altitude</b>
1	32	0.4	10	30
2	32	0.4	15	35
3	32	0.4	20	40
4	32	0.4	25	50
5	32	0.4	30	80

You could then plot the data on a graph with the pressure on the horizontal axis and the altitude on the vertical axis. The next series of experiments might increase or decrease the amount of water used and then repeat the varying pressures.

Other interesting experiments might hold the pressure constant and vary the amount of water used. Try changing the rocket's length, or the rocket's weight; the combinations are almost endless! Just remember, keep records, plot graphs, and only change one parameter at time!

**Extensions:**

You may want to have competitions with friends to see whose rocket will go the highest, stay in the air the longest, carry the most weight, or land closest to a specified spot.

**Time Allocation:**

This product arrives ready to use. Allow 10 minutes to set up and launch the first soda bottle. Individual experiment times will vary depending on the complexity of the rocket design and the methods of instruction, but ordinarily no single activity will exceed one class period.

19 18 17 16 15 14 13 12 11 10 9 8 7 6 5 4 3 2 1




1 - Sight rocket along top edge of altimeter.  
2 - When rocket is at peak altitude, pinch the plumb line against the graph paper.  
3 - look along the plumb line for a convenient intersection.  
4 - Determine X and Y from this intersection.

**Height = Base \* (Y / X)**

Glue this sheet onto a backing of cardboard.

Punch hole for string here



X  
1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15

# -Rocket Altimeter-