

Health Occupations



LF01095H

Life/form® Blood Pressure Simulator

This Blood Pressure Simulator helps resolve the uncertainties common in teaching students to take blood pressure. This lifelike simulator allows the presetting of values for both systolic and diastolic pressures. Provides an excellent means to practice listening to and distinguishing blood pressure sounds prior to actual clinical experience. It is possible to audibly discern the five Korotkoff phases. Many times when working with a live subject, pressures are difficult to auscultate, making accurate evaluation of student proficiency almost impossible, and undermining the student's self-confidence. With this realistic unit, the student can find the preset results and the instructor can unfailingly know if the student has performed the procedure accurately. Teacher can easily calibrate unit for use with any sphygmomanometer. The Blood Pressure Simulator includes the Life/form® Arm Replica with embedded speaker, sphygmomanometer, electronic control unit (six "AA" batteries included), instruction booklet, and hard carry case. The simulator does not include a stethoscope. For auscultatory blood pressure only. 23" L x 4" W x 4" D. Five-year warranty. Sh. wt. 10 lbs.

LF01095H — \$965.00

ISO 9001 & ISO 13485



SB40247H

SB27332H

NEW Aneroid Fanny Pack

Includes color-matched Sprague-Rappaport type stethoscope with accessory kit, no-stop manometer gauge, color-matched nylon cuff latex inflation system, stainless steel bandage scissors, reusable pen light, tape measure, and stethoscope name tag. Ten-year manometer warranty. Sh. wt. 2.50 lbs.

SB40246H Burgundy SB40248H Royal Blue SB40250H Black
SB40247H Purple SB40249H Teal
Each — \$60.00

Teaching Sphygmomanometer Set

Designed as an instructional set. Includes a dual-head teaching stethoscope plus a popular aneroid sphygmomanometer. To be used for simultaneous listening by teacher and student or two students. Excellent device for assisting students in the detection of systolic and diastolic heart sounds. Sh. wt. 0.81 lbs.

SB27332H — \$30.50

NEW Fingertip Pulse Oximeter

Safely monitor oxygen levels and pulse rate. Accurately measure oxygen saturation of arterial hemoglobin and pulse rate. Accommodates finger sizes from 1/4" to 1 1/16" dia., or 33-lb. to 242-lb. users. Large LCD display with red digital readout displays blood-oxygen saturation (SpO2), pulse rate, and pulse bar. One-button operation, low battery indicator, and automatic shut-off. Includes lanyard and detailed English/Spanish instructions. Requires two "AAA" batteries (included). Latex free. One-year limited warranty. Sh. wt. 0.31 lbs.

SB48022H — \$29.95



NEW Nasco's Vision Challenges Simulation Glasses Set

Create the realities of six common age-related visual impairments: cataracts, glaucoma, macular degeneration, retinopathy, hemianopsia, and detached retina. Designed to be easily worn over glasses. Use glasses to attempt common activities of daily living (filling out forms, reading the newspaper, mixing cooking ingredients, etc.). Includes six pairs of glasses and user guide with visual challenges worksheets. Sh. wt. 0.19 lbs.

SB50159H — \$29.50



FEATURED PRODUCT



5-Minute Health Sciences Activities

From Learning ZoneXpress™ Gr. 7+ Fifty quick and effective activities that get students thinking about all aspects of health care and sciences. Filled with creative, fun, and thought-provoking activity ideas, including medical ethics, health care services and careers, medical equipment, body systems, nutrition, and more. 50 pages, 8 1/2" x 5 1/2". Spiral bound.

©2014. Sh. wt. 0.88 lbs.
WA32531H — \$29.95 5+ \$27.95

DVD



Health Science Careers

Gr. 7-12 This program explores a variety of health science career options, with varying levels of education required, including pharmacy technician, home health aid, pediatrician, medical administrative support, radiographer, physical therapist and physical therapist assistant, and certified nursing assistant. Teaching materials included. 20 minutes. ©2010. Sh. wt. 0.31 lbs.

WA28718H DVD — \$79.95

Software


Anatomy of Health Careers

Created by Masters Prepared Educators. Second Edition.

Zoom in and enjoy this fabulous CD-ROM! An incredible compilation that contains 62 of the most in demand health careers. Divided into 15 "user-friendly" sections based on occupational characteristics. Each career page contains an informative interview with a professional in each of the occupations, typical places of employment, entry-level or average starting salary information, requirements to enter the careers, a brief description on the career and list of specializations, extensive financial aid information, professional organizations to contact for further information, and links to websites. Sh. wt. 0.31 lbs.

SB40265H CD-ROM — \$99.00

SYSTEM REQUIREMENTS: Windows® 98 or higher, 8 MB RAM; Macintosh® OS 9.1 and OS X®, Adobe® Acrobat® Reader.



Quick & Easy Medical Terminology

By Peggy C. Leonard, BA, MT, MEd.

Eighth Edition.

Read, write, and speak medical terms in no time. Ideal for one-semester courses and self-study, this clear, concise text covers disorders, diagnosis, and treatment with a strong focus on anatomy for an integrated discussion on each body system. Small bits of information are presented then immediately tested for comprehension using a combination of programmed and traditional review questions. With new case studies, terms, and photos, plus help from the companion Evolve site, Quick & Easy Medical Terminology provides you with simple tools for good communication in the health care environment. Now includes a brief explanation and focus on basic word-building skills, moderate anatomy and physiology coverage, programmed learning, review exercises integrated throughout lesson material, medical reports with review exercises, quick tips throughout the text, and Student Evolve resources. 456 pages, 7 1/2" x 9 1/4". Spiral-bound. ©2016. Sh. wt. 2 lbs.

SB47365H Book & CD-ROM — \$63.95

SYSTEM REQUIREMENTS: Windows® XP(SP2)/Vista, Intel Pentium® IV processor, 512 MB RAM; Macintosh® OS 10.4 and 10.5, Intel-based processor, 512 MB RAM.