

8165

PE09718
PE09719
V# 022530



SAFETY DATA SHEET

Issuing Date 18-Jun-2013

Revision Date 18-Jun-2013

Revision Number 0

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

GHS product identifier

Product Name Firm Grip, QDA, Tuf-Skin

Other means of identification

Product Code(s) 061027, 171531, 204027, 204028, 204033, 203531

UN-Number UN 1950

Synonyms None

Recommended use of the chemical and restrictions on use

Recommended Use Apply to skin prior to taping, or to fingertips for improved grip

Uses advised against No information available

Supplier's details

Supplier Address
Cramer Products, Inc.
153 W. Warren
Gardner, KS 66030
TEL: 913.856.7511

Emergency telephone number

Emergency Telephone Number 913 856.7511

2. HAZARDS IDENTIFICATION


Classification

This chemical is considered hazardous according to the OSHA Hazard Communication Standard 2012 (29 CFR 1910.1200)

Serious Eye Damage/Eye Irritation	Category 2
Skin Sensitization	Category 1
Specific Target Organ Systemic Toxicity (Single Exposure)	Category 3
Simple Asphyxiant	Yes
Flammable aerosols	Category 2

GHS Label elements, including precautionary statements

Emergency Overview

Signal Word	Warning
Hazard Statements	
<ul style="list-style-type: none"> • Causes serious eye irritation • May cause an allergic skin reaction • May cause drowsiness or dizziness • May displace oxygen and cause rapid suffocation • Flammable aerosol 	
	
Appearance Clear to amber	Physical State Aerosol.
	Odor IPA, rosin

Precautionary Statements**Prevention**

- Keep away from heat/sparks/open flames/hot surfaces - No smoking
- Pressurized container: Do not pierce or burn, even after use
- Do not spray on an open flame or other ignition source
- Wear protective gloves
- Wear eye/face protection
- Avoid breathing dust/fume/gas/mist/vapors/spray
- Wash face, hands and any exposed skin thoroughly after handling
- Contaminated work clothing should not be allowed out of the workplace
- Use only outdoors or in a well-ventilated area

General Advice

- None

Eyes

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
- If eye irritation persists: Get medical advice/attention.

Skin

- IF ON SKIN: Wash with plenty of soap and water
- If skin irritation or rash occurs: Get medical advice/attention
- Wash contaminated clothing before reuse.

Inhalation

- IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
- Call a POISON CENTER or doctor/physician if you feel unwell

Storage

- Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F
- Store in a well-ventilated place. Keep container tightly closed
- Store locked up.

Disposal

- Dispose of in accordance with local/regional/national regulations.

Hazard Not Otherwise Classified (HNOC)

May cause frostbite.

Other information

No information available.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	Weight %	Trade secret
Isopropyl alcohol	67-63-0	30-60	*
Petroleum gases, liquified, sweetened	68476-86-8	15-40	*
Resin	8050-09-7	10-30	*
Benzoin Tincture IPA	-	1-5	*

*The exact percentage (concentration) of composition has been withheld as a trade secret.

4. FIRST AID MEASURES

Description of necessary first-aid measures

- Eye Contact** Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Keep eye wide open while rinsing. Consult a physician.
- Skin Contact** Wash off immediately with soap and plenty of water removing all contaminated clothes and shoes. In case of contact with liquefied gas, thaw frosted parts with lukewarm water. If symptoms persist, call a physician.
- Inhalation** Move to fresh air. Administer oxygen if breathing is difficult and you are trained. If breathing is irregular or stopped, administer artificial respiration. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Call a physician or Poison Control Center immediately.
- Ingestion** Do NOT induce vomiting. Rinse mouth. Drink plenty of water. Never give anything by mouth to an unconscious person. Consult a physician.

Most important symptoms/effects, acute and delayed

Most important Symptoms/Effects Difficulty in breathing. Dizziness. Drowsiness. Itching. Rashes. Frostbite. Red eyes, Tearing.

Indication of immediate medical attention and special treatment needed, if necessary

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media CAUTION: Use of water spray when fighting fire may be inefficient.

Specific Hazards Arising from the Chemical

Vapors may form explosive mixtures with air. Vapors may travel to source of ignition and flash back. Thermal decomposition can lead to release of irritating and toxic gases and vapors. In the event of fire and/or explosion do not breathe fumes.

Explosion Data

Sensitivity to Mechanical Impact	Yes.
Sensitivity to Static Discharge	Yes.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions Remove all sources of ignition. Take precautionary measures against static discharges. All equipment used when handling the product must be grounded. Evacuate personnel to safe areas. Keep people away from and upwind of spill/leak. Ensure adequate ventilation. Use personal protective equipment. Do not breathe vapors or spray mist. Do not get in eyes, on skin, or on clothing. Contents under pressure. Do not puncture or incinerate cans.

Environmental Precautions

Environmental Precautions See Section 12 for additional Ecological Information

Methods and materials for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Soak up with inert absorbent material. Take up mechanically and collect in suitable container for disposal. Ground and bond containers when transferring material.

7. HANDLING AND STORAGE

Precautions for safe handling

Handling Contents under pressure. Do not puncture or incinerate cans. Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharges. Wear personal protective equipment. Ensure adequate ventilation. Do not breathe vapors or spray mist. Do not get in eyes, on skin, or on clothing.

Conditions for safe storage, including any incompatibilities

Storage Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from open flames, hot surfaces and sources of ignition. Keep away from direct sunlight.

Incompatible Products Strong oxidizing agents. Acids. Acid anhydrides. Aluminum. Halogenated compounds.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Control parameters

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Isopropyl alcohol 67-63-0	STEL: 400 ppm TWA: 200 ppm	TWA: 400 ppm TWA: 980 mg/m ³ (vacated) TWA: 400 ppm (vacated) TWA: 980 mg/m ³ (vacated) STEL: 500 ppm (vacated) STEL: 1225 mg/m ³	IDLH: 2000 ppm 15% LEL TWA: 980 mg/m ³ TWA: 400 ppm STEL: 500 ppm STEL: 1225 mg/m ³
Rosin 8050-09-7		(vacated) TWA: 0.1 mg/m ³ Formaldehyde	TWA: 0.1 mg/m ³ Formaldehyde

Appropriate engineering controls

Engineering Measures
 Showers
 Eyewash stations
 Ventilation systems

Individual protection measures, such as personal protective equipment

Eye/Face Protection Tightly fitting safety goggles. If splashes are likely to occur, wear: Face-shield.
Skin and Body Protection Impervious gloves. Impervious clothing. Boots.
Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or smoke when using this product. Keep away from food, drink and animal feeding stuffs. Wash hands before breaks and immediately after handling the product. Remove and wash contaminated clothing before re-use. Provide regular cleaning of equipment, work area and clothing.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State Odor	Aerosol IPA, rosin	Appearance Odor Threshold	Clear to amber No information available
Property	Values	Remarks/ - Method	
pH	No data available	None known	
Melting Point/Range	No data available	None known	
Boiling Point/Boiling Range	No data available	None known	
Flash Point	No data available	None known	
Evaporation rate	No data available	None known	
Flammability (solid, gas)	No data available	None known	
Flammability Limits in Air upper flammability limit	No data available		
lower flammability limit	No data available		
Vapor Pressure	No data available	None known	
Vapor Density	No data available	None known	
Relative Density	No data available	None known	
Specific Gravity	No data available	None known	
Water Solubility	Insoluble in water.	None known	
Solubility in other solvents	No data available	None known	
Partition coefficient: n-octanol/water	No data available	None known	
Autoignition Temperature	No data available	None known	
Decomposition Temperature	No data available	None known	
Viscosity	No data available	None known	
Flammable Properties	Flammable aerosol.		
Explosive Properties	No data available		
Oxidizing Properties	No data available		
Other information			
VOC Content (%)	No data available		

10. STABILITY AND REACTIVITY

Reactivity

No data available.

#445 P.009/096
 11/27/2017 17:15
 913 884+5626
 From: CRAMER SPORTS MED

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

None under normal processing.

Conditions to avoid

Keep away from open flames, hot surfaces and sources of ignition. Keep away from direct sunlight.

Incompatible materials

Strong oxidizing agents, Acids, Acid anhydrides, Aluminium, Halogenated compounds.

Hazardous decomposition products

Carbon monoxide (CO), Carbon dioxide (CO₂).

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information	No acute toxicity information is available for this product. Description of possible hazardous to health effects is based on experience and/or toxicological characteristics of several ingredients.
Inhalation	May cause drowsiness and dizziness. Deliberately inhaling or exposure to high concentrations may cause simple asphyxiation (suffocation by displacing the oxygen in the air).
Eye Contact	Irritating to eyes.
Skin Contact	May cause sensitization by skin contact. Contact with product may cause frostbite.
Ingestion	Not an expected route of exposure. May be harmful if swallowed.

Toxicology data for the components

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Isopropyl alcohol	= 4396 mg/kg (Rat)	12800 mg/kg (Rat) 12670 mg/kg (Rabbit)	72.6 mg/L (Rat) 4 h
Rosin	= 3 mg/kg (Rat)	> 2500 mg/kg (Rabbit)	-

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms	Irritation, Eye and skin redness, tearing, hives, blurry vision. Vapors may cause drowsiness and dizziness. Inhalation of high vapor concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.
-----------------	--

Delayed and immediate effects and also chronic effects from short and long term exposure

Sensitization	May cause sensitization by skin contact.
Mutagenic Effects	No information available.
Carcinogenicity	The table below indicates whether each agency has listed any ingredient as a carcinogen.

Chemical Name	ACGIH	IARC	NTP	OSHA
Isopropyl alcohol		Group 3		X

IARC: (International Agency for Research on Cancer)
 Group 3: Not Classifiable as to its Carcinogenicity to Humans
OSHA: (Occupational Safety & Health Administration)
 X - Present

Reproductive Toxicity	No information available.
STOT - single exposure	May cause drowsiness or dizziness.

#445 P.010/096

11/27/2017 17:16

913 884+5626

From: CRAMER SPORTS MED

STOT - repeated exposure No information available.
Aspiration Hazard No information available.

Numerical measures of toxicity - Product
Acute Toxicity No information available.

12. ECOLOGICAL INFORMATION

Ecotoxicity

The environmental impact of this product has not been fully investigated

Chemical Name	Toxicity to Algae	Toxicity to Fish	Toxicity to Microorganisms	Daphnia Magna (Water Flea)
Isopropyl alcohol	EC50 96 h: > 1000 mg/L (Desmodesmus subspicatus) EC50 72 h: > 1000 mg/L (Desmodesmus subspicatus)	LC50 96 h: = 11130 mg/L static (Pimephales promelas) LC50 96 h: = 9640 mg/l flow-through (Pimephales promelas) LC50 96 h: > 1400000 µg/L (Lepomis macrochirus)		EC50 48 h: = 13299 mg/L (Daphnia magna)
Resin	EC50 72 h: = 400 mg/L (Desmodesmus subspicatus)		EC50 = 31.5 mg/L 30 min	EC50 48 h: 3.8 - 5.4 mg/L (Daphnia magna)

Persistence and Degradability No information available.

Bioaccumulation No information available.

Chemical Name	Log Pow
Isopropyl alcohol	0.05
Petroleum gases, liquified, sweetened	2.8

Other Adverse Effects
No information available

13. DISPOSAL CONSIDERATIONS

Waste Disposal Methods This material, as supplied, is a hazardous waste according to federal regulations (40 CFR 261).

Contaminated Packaging Do not re-use empty containers.

US EPA Waste Number D001

This product contains one or more substances that are listed with the State of California as a hazardous waste.

Chemical Name	California Hazardous Waste
Isopropyl alcohol	Toxic Ignitable

14. TRANSPORT INFORMATION

Note: Limited quantity

DOT
UN-Number UN 1950
Proper shipping name Aerosols, flammable, n.o.s.
Hazard Class 2.1

#445 P.011/096

11/27/2017 17:17

913 884+5626

From: CRAMER SPORTS MED

15. REGULATORY INFORMATION

International Inventories

TSCA Exempt
DSL Exempt

Legend

TSCA - United States Toxic Substances Control Act; Section 8(h) Inventory
DSL/NDL - Canadian Domestic Substances List/Non-Domestic Substances List

U.S. Federal Regulations

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Table with 4 columns: Chemical Name, CAS-No, Weight %, SARA 313 - Threshold Values %. Row 1: Isopropyl alcohol, 67-83-0, 54.15, 1.0

SARA 311/312 Hazard Categories

Acute Health Hazard Yes
Chronic Health Hazard Yes
Fire Hazard Yes
Sudden Release of Pressure Hazard Yes
Reactive Hazard No

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

U.S. State Right-to-Know Regulations

Table with 6 columns: Chemical Name, New Jersey, Massachusetts, Pennsylvania, Illinois, Rhode Island. Rows: Isopropyl alcohol, Rosin

U.S. EPA Label Information

EPA Pesticide Registration Number Not applicable

16. OTHER INFORMATION

NFPA Health Hazard 2 Flammability 4 Instability 0 Physical and Chemical Hazards -
HMIS Health Hazard 2 Flammability 4 Physical Hazard 0 Personal Protection X

#445 P.012/096

11/27/2017 17:17

913 884+5626

From: CRAMER SPORTS MED

WPS-CPI-002 - Firm Grip, ODA, Tuf-Skin

Revision Date 18-Jun-2013

Prepared By	Product Stewardship 23 British American Blvd. Latham, NY 12110 1-800-572-6501
Issuing Date	18-Jun-2013
Revision Date	18-Jun-2013
Revision Note	Initial Release.

General Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of Safety Data Sheet