




## BLOOD CLOTS, PERINEAL, AND BLOOD POOL

### **Cleaning Instructions**

- 1.** To ensure life of product, wash with warm soapy water and completely air dry before returning to bag.
- 2.** Clot material may stain manikin; contact with ink from printed material and ballpoint pen will also stain manikin.
- 3.** Clots are to be used in simulation only.



**Small Blood Clot**  
**LF00710U**



**Medium Blood Clot**  
**LF00709U**



**Large Blood Clot**  
**LF00711U**



**Perineal Hemorrhage**  
**LF00712U**



**Hemorrhage Blood Pool**  
**LF00713U**

**Complete Set of Clots and Hemorrhages**  
**LF00714U**

ISO 9001 and  
ISO 13485

*Nasco* **Plastics**  
ISO 9001 and ISO 13485

Products are manufactured using a quality management system certified as being in conformity with ISO 9001 and ISO 13485.