

**Date:** April 1, 2014

**From:** Victor Lefaive – Nasco Plastics


**Subject:** Product Improvement – Bariatric CPR Manikin  
(Fat Old Fred) – LF03749 and LF03750



Nasco is pleased to announce the improved redesign of internal structures of the **Bariatric CPR Manikin LF03749 and LF03750** (Manikin formerly known as “Fat Old Fred”).

We have upgraded the foam sternum insert (LF03777U) to have a PVC skin to protect it and prolong the life of the product. The new insert has been tested, *without failure*, to 24,000 compressions – a *significant* improvement to the original design. All new Bariatric CPR Manikins will have this feature and all existing stock has been inspected and upgraded to the new version.

**LF03777** – Foam Sternum Insert  
Replacement for LF03749 & LF03750  
(Old version without skin)



**LF03777** – Foam Sternum  
Insert Replacement for  
LF03749 & LF03750  
(New PVC skin to prolong  
product life)



If you have any further questions, please contact:

**Nasco - Fort Atkinson, WI**

Phone: 1-800-558-9595 | 1-920-563-2446

Fax: 1-920-563-8296

E-mail: [info@eNasco.com](mailto:info@eNasco.com) or  
[export@eNasco.com](mailto:export@eNasco.com)

Weekdays: 7 a.m. - 6 p.m.

Saturday: 8 a.m. - 12 noon

**Nasco - Modesto, CA**

Phone: 1-800-558-9595 | 1-209-545-1600

Fax: 1-209-545-1669

E-mail: [modestocs@eNasco.com](mailto:modestocs@eNasco.com)

Weekdays: 7 a.m. - 6 p.m.

Saturday: 8 a.m. - 3 p.m.

Thank you,

Victor Lefaive

Quality Assurance Coordinator

Nasco Plastics

New style foam sternum  
insert, quality tested to  
**24,000 compressions**  
*without failure*, prolongs  
product life!

