

To: Sales and Customer Service
From: Victor Lefaive
Subject: Catheter Obstruction Tool and User Assistance Videos
Date: 8//2014

Nasco Plastics has updated the **LF00855U Male Catheterization Simulator** and the **LF00856U Female Catheterization Simulator** to include an obstruction tool that removes blockages caused by lubricant buildup. Customers who already own catheter simulators may purchase **LF00844 Catheter Obstruction Tool Kit** separately as needed.



Tool used to clear obstructions caused by built-up/dried lubricant.

Nasco has also created user help **videos** to assist customers with the proper use the catheterization simulators as well as the obstruction tool.

The videos can be found on the product page for the male and female catheterization simulators (LF00855 & LF00856) on the Nasco website – www.eNasco.com/healthcare/.



Click the images above to access the product page, videos, and newly revised instruction manuals.

Thank you,

Victor Lefaive
Quality Assurance Coordinator
Nasco Plastics