


Date: 11/20/13

To: Domestic Distributors

From: Victor Lefaive – Nasco Plastics

Subject: Product Improvement – Infant *CRiSis*TM tubing update

Nasco Plastics has updated the vein tubing for the IV arms and legs on the infant *CRiSis*TM manikins to improve product performance. The **shrink wrap** previously used to encase the vein tubing has been eliminated to increase the flow of fluid through the simulator.



Additionally, the **tubing clamps** will no longer be attached to the veins when packaged for shipping; rather they will simply accompany the product in order to prevent kinking of the veins during shipment/storage.

Customers should notice a dramatic improvement in the functionality of the product. All stock has been pulled from the Nasco warehouse and upgraded.

Affected Products:

[LF03635U](#) Complete PALS Add-on Kit

[LF03637U](#) Infant IV Arm

[LF03640U](#) Infant IV Leg Replacement Veins

LF03645U PALS IV Arm, Leg
and Umbilicus Chest Skin

[LF03636U](#) Infant IV Leg

[LF03639U](#) IV Leg Replacement Skin/Veins

[LF03641U](#) IV Arm Replacement Skin/Veins

[LF03709U](#) Complete Infant *CRiSis*TM Manikin

[LF03700U](#) Intermediate Infant *CRiSis*TM Manikin

If you have any further questions, please contact:

Nasco - Fort Atkinson, WI

Phone: 1-800-558-9595 | 1-920-563-2446

Fax: 1-920-563-8296

E-mail: info@enasco.com

Weekdays: 7 a.m. - 6 p.m.

Saturday: 8 a.m. - 12 noon

Nasco - Modesto, CA

Phone: 1-800-558-9595 | 1-209-545-1600

Fax: 1-209-545-1669

E-mail: modesto@enasco.com

Weekdays: 7 a.m. - 6 p.m.

Saturday: 8 a.m. - 3 p.m.

Thank you,

Victor Lefaive, Quality Assurance Coordinator
Nasco Plastics